


Anesthetic Products




Standard Kits include:

1. Central Venous Catheter
2. Guide-wire
3. Vessel Dilator
4. Clamp
5. Fastener: Catheter Clamp
6. Introducer Needle
7. Introducer Syringe
8. Injection Needle
9. Injection Cap

Optional Compound Kits include:

1. Central Venous Catheter
2. Standard Kit Accessories
3. 5ml Syringe
4. Surgical Gloves
5. Surgical Pledget
6. Surgery Sheet
7. Surgery Towel
8. Sterile Brush
9. Gauze Pad
10. Suture of Needle
11. Wound Dressing
12. Scalpel

Central Venous Catheter



Movable clamp allows anchorage at puncture site regardless of depth of catheter, which minimizes trauma and irritation to the puncturing site. Depth marking assists in accurate placement of central venous catheter from the right or left subclavian or jugular vein. Soft tip reduces trauma to vessel, minimizing vessel erosion, hemothorax and cardiac tamponade. Single, double, triple and quad lumen is available for choice. Radiopacity facilitates confirmation of catheter placement. The multi-lumen version's tips are more radiopaque, which makes it easy to confirm the placement of fluoroscopic tip. Vessel dilator ensures "super soft" catheters to be easily placed percutaneously.

Optional




Triangle Hub

Alternative



Straight or Y-shaped Type Puncture Needle.



Nitinol or Stainless Steel Guide-wire with Flexible J-tip



Sharp Scalpel




ECG Cable

Compound-packaging Tray



Various Accessories for Choice

Order Information

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV-1221	Single	22Ga	5	4F×5cm	0.018"×50cm	Straight/Y Needle 20Ga 
FV-1222	Single	22Ga	8	4F×5cm	0.018"×50cm	
FV-1223	Single	22Ga	10	4F×5cm	0.018"×50cm	
FV-1224	Single	22Ga	13	4F×5cm	0.018"×50cm	
FV-122915	Single	22Ga	15	4F×5cm	0.018"×50cm	
FV-1225	Single	22Ga	16	4F×5cm	0.018"×50cm	
FV-1226	Single	22Ga	20	4F×5cm	0.018"×60cm	
FV-1228	Single	22Ga	30	4F×5cm	0.018"×70cm	
FV-1321	Single	20Ga	5	4F×5cm	0.021"×50cm	Straight/Y Needle 20Ga 
FV-1322	Single	20Ga	8	4F×5cm	0.021"×50cm	
FV-1323	Single	20Ga	10	4F×5cm	0.021"×50cm	
FV-1324	Single	20Ga	13	4F×5cm	0.021"×50cm	
FV-132915	Single	20Ga	15	4F×5cm	0.021"×50cm	
FV-1325	Single	20Ga	16	4F×5cm	0.021"×50cm	
FV-1326	Single	20Ga	20	4F×5cm	0.021"×60cm	
FV-1328	Single	20Ga	30	4F×5cm	0.021"×70cm	
FV-1421	Single	18Ga	5	5.5F×6cm	0.025"×50cm	Straight/Y Needle 19Ga 
FV-1422	Single	18Ga	8	5.5F×6cm	0.025"×50cm	
FV-1423	Single	18Ga	10	5.5F×6cm	0.025"×50cm	
FV-1424	Single	18Ga	13	5.5F×6cm	0.025"×50cm	
FV-142915	Single	18Ga	15	5.5F×6cm	0.025"×50cm	
FV-1425	Single	18Ga	16	5.5F×6cm	0.025"×50cm	
FV-1426	Single	18Ga	20	5.5F×6cm	0.025"×60cm	
FV-1428	Single	18Ga	30	5.5F×6cm	0.025"×70cm	
FV-1521	Single	16Ga	5	6.5F×10cm	0.032"×50cm	Straight/Y Needle 18Ga 
FV-1522	Single	16Ga	8	6.5F×10cm	0.032"×50cm	
FV-1523	Single	16Ga	10	6.5F×10cm	0.032"×50cm	
FV-1524	Single	16Ga	13	6.5F×10cm	0.032"×50cm	
FV-152915	Single	16Ga	15	6.5F×10cm	0.032"×50cm	
FV-1525	Single	16Ga	16	6.5F×10cm	0.032"×50cm	
FV-1526	Single	16Ga	20	6.5F×10cm	0.032"×60cm	
FV-1528	Single	16Ga	30	6.5F×10cm	0.032"×70cm	

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV-1624	Single	14Ga	13	8.5F×10cm	0.032"×50cm	Straight/Y Needle 18Ga 
FV-162915	Single	14Ga	15	8.5F×10cm	0.032"×50cm	
FV-1625	Single	14Ga	16	8.5F×10cm	0.032"×50cm	
FV-1626	Single	14Ga	20	8.5F×10cm	0.032"×60cm	
FV-1628	Single	14Ga	30	8.5F×10cm	0.032"×70cm	Straight/Y Needle 20Ga 
FV/FC-2421	Double	4Fr	5	5.5F×6cm	0.018"×50cm	
FV/FC-2422	Double	4Fr	8	5.5F×6cm	0.018"×50cm	
FV/FC-2423	Double	4Fr	10	5.5F×6cm	0.018"×50cm	
FV/FC-2424	Double	4Fr	13	5.5F×6cm	0.018"×50cm	Straight/Y Needle 20Ga 
FV/FC-242915	Double	4Fr	15	5.5F×6cm	0.018"×50cm	
FV/FC-2425	Double	4Fr	16	5.5F×6cm	0.018"×50cm	
FV/FC-2426	Double	4Fr	20	5.5F×6cm	0.018"×60cm	
FV/FC-2428	Double	4Fr	30	5.5F×6cm	0.018"×70cm	Straight/Y Needle 20Ga 
FV/FC-2521	Double	5Fr	5	6F×7cm	0.021"×50cm	
FV/FC-2522	Double	5Fr	8	6F×7cm	0.021"×50cm	
FV/FC-2523	Double	5Fr	10	6F×7cm	0.021"×50cm	
FV/FC-2524	Double	5Fr	13	6F×7cm	0.021"×50cm	Straight/Y Needle 18Ga 
FV/FC-252915	Double	5Fr	15	6F×7cm	0.021"×50cm	
FV/FC-2525	Double	5Fr	16	6F×7cm	0.021"×50cm	
FV/FC-2526	Double	5Fr	20	6F×7cm	0.021"×60cm	
FV/FC-2528	Double	5Fr	30	6F×7cm	0.021"×70cm	Straight/Y Needle 18Ga 
FV/FC-2724	Double	7Fr	13	8.5F×10cm	0.035"×50cm	
FV/FC-272915	Double	7Fr	15	8.5F×10cm	0.035"×50cm	
FV/FC-2725	Double	7Fr	16	8.5F×10cm	0.035"×50cm	
FV/FC-2726	Double	7Fr	20	8.5F×10cm	0.035"×60cm	Straight/Y Needle 18Ga 
FV/FC-2728	Double	7Fr	30	8.5F×10cm	0.035"×70cm	
FV-2824	Double	8Fr	13	9F×10cm	0.035"×50cm	
FV-282915	Double	8Fr	15	9F×10cm	0.035"×50cm	
FV-2825	Double	8Fr	16	9F×10cm	0.035"×50cm	Straight/Y Needle 18Ga 
FV-2826	Double	8Fr	20	9F×10cm	0.035"×50cm	
FV-2828	Double	8Fr	30	9F×10cm	0.035"×50cm	

Order Information

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV-2921	Double	8.5Fr	5	9F×10cm	0.035"×50cm	Straight/Y Needle 18Ga 
FV-2922	Double	8.5Fr	8	9F×10cm	0.035"×50cm	
FV-2923	Double	8.5Fr	10	9F×10cm	0.035"×50cm	
FV-2924	Double	8.5Fr	13	9F×10cm	0.035"×50cm	
FV-2925	Double	8.5Fr	16	9F×10cm	0.035"×50cm	
FV-2926	Double	8.5Fr	20	9F×10cm	0.035"×60cm	
FV-2927	Double	8.5Fr	25	9F×10cm	0.035"×60cm	
FV-2928	Double	8.5Fr	30	9F×10cm	0.035"×70cm	
FV/FC-3421	Triple	4.5Fr	5	6F×5cm/6F×7cm	0.021"×50cm	Straight/Y Needle 20Ga 
FV/FC-3422	Triple	4.5Fr	8	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3423	Triple	4.5Fr	10	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3424	Triple	4.5Fr	13	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-342915	Triple	4.5Fr	15	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3425	Triple	4.5Fr	16	6F×5cm/6F×7cm	0.021"×50cm	
FV/FC-3426	Triple	4.5Fr	20	6F×5cm/6F×7cm	0.021"×60cm	
FV/FC-3428	Triple	4.5Fr	30	6F×5cm/6F×7cm	0.021"×70cm	
FV/FC-3521	Triple	5.5Fr	5	6.5F×10cm	0.021"×50cm	Straight/Y Needle 20Ga 
FV/FC-3522	Triple	5.5Fr	8	6.5F×10cm	0.021"×50cm	
FV/FC-3523	Triple	5.5Fr	10	6.5F×10cm	0.021"×50cm	
FV/FC-3524	Triple	5.5Fr	13	6.5F×10cm	0.021"×50cm	
FV/FC-352915	Triple	5.5Fr	15	6.5F×10cm	0.021"×50cm	
FV/FC-3525	Triple	5.5Fr	16	6.5F×10cm	0.021"×50cm	
FV/FC-3526	Triple	5.5Fr	20	6.5F×10cm	0.021"×60cm	
FV/FC-3528	Triple	5.5Fr	30	6.5F×10cm	0.021"×70cm	
FV/FC-3724	Triple	7Fr	13	8.5F×10cm	0.035"×50cm	Straight/Y Needle 18Ga 
FV/FC-372915	Triple	7Fr	15	8.5F×10cm	0.035"×50cm	
FV/FC-3725	Triple	7Fr	16	8.5F×10cm	0.035"×50cm	
FV/FC-3726	Triple	7Fr	20	8.5F×10cm	0.035"×60cm	
FV/FC-3728	Triple	7Fr	30	8.5F×10cm	0.035"×70cm	Straight/Y Needle 18Ga 

Product Code	Lumen	Catheter Size	Length (cm)	Dilator	Guidewire	Introducer Needle
FV/FC-3924	Triple	8.5Fr	13	9F×10cm	0.035"×50cm	Straight/Y Needle 18Ga 
FV/FC-392915	Triple	8.5Fr	15	9F×10cm	0.035"×50cm	
FV/FC-3925	Triple	8.5Fr	16	9F×10cm	0.035"×50cm	
FV/FC-3926	Triple	8.5Fr	20	9F×10cm	0.035"×60cm	
FV/FC-3928	Triple	8.5Fr	30	9F×10cm	0.035"×70cm	Straight/Y Needle 18Ga 
FV-4921	Quad	8.5Fr	5	9F×10cm	0.035"×50cm	
FV-4922	Quad	8.5Fr	8	9F×10cm	0.035"×50cm	
FV-4923	Quad	8.5Fr	10	9F×10cm	0.035"×50cm	
FV-4924	Quad	8.5Fr	13	9F×10cm	0.035"×50cm	
FV-492915	Quad	8.5Fr	15	9F×10cm	0.035"×50cm	
FV-4925	Quad	8.5Fr	16	9F×10cm	0.035"×50cm	Straight/Y Needle 18Ga 
FV-4926	Quad	8.5Fr	20	9F×10cm	0.035"×50cm	
FV-4928	Quad	8.5Fr	30	9F×10cm	0.035"×50cm	

Remarks:

1. The product codes in this chart are for CVC kits with straight introducer needle configuration.
2. CVC kits with Y-shaped introducer needle configuration are also available in adult types. The product code of CVC kit with Y-shaped introducer needle configuration is to add letter "Y" after the product code of CVC kit with straight introducer needle configuration, for example, "FC-2726Y". Before you order, please clarify which configuration you request, either straight or Y-shape introducer needle.
3. Pediatric types of CVC include single lumen 22Ga, 20Ga, 18Ga ; double lumen 4Fr, 5Fr ; triple lumen 4.5Fr and triple lumen 5.5Fr.
4. Adult types of CVC include single lumen 16Ga, 14Ga; double lumen 7Fr, 8Fr, 8.5Fr; triple lumen 7Fr, 8.5Fr and quad lumen 8.5Fr.
5. "FV" in the product code for single lumen CVC means with cylinder hub.
6. "FC" in the product code for double lumen, triple lumen, quad lumen CVC means with triangular hub. "FV" in the product code for double lumen, triple lumen, quad lumen CVC means with round hub.