

# MASKS



## ANESTHESIA MASKS

- Designed for anesthesia, respiratory or resuscitation
- Leak-proof, soft and comfortable seal cushion in facial contoured design
- Color coded hook ring for easy size identification
- Single patient use to prevent cross contamination
- Light weight for patient comfort
- Standard and scented types
- Latex-free
- Available in different sizes
- Meets ISO Standards
- Transparent to allow easy observation



## Economy Masks

	<b>Code</b>	<b>Size</b>	<b>Weight</b>	<b>Dead Space</b>	<b>Qty/Box</b>
Respiratory, anesthesia mask #5	8050	Large	49,4 gr	255 ml	25
Respiratory, anesthesia mask #4	8040	Medium	38,6 gr	200 ml	25
Respiratory, anesthesia mask #3	8030	Small	28,9 gr	140 ml	25
Respiratory, anesthesia mask #2	8020	Child	23,6 gr	92 ml	50
Respiratory, anesthesia mask #1	8010	Infant	13,3 gr	51 ml	50
Respiratory, anesthesia mask #0	8000	Neo-nate	9,8 gr	28 ml	50