

STYLETS

FLEXISLIP INTUBATION STYLET

flexible metal stylet with plastic coating, smooth surface, curved handle

- for tracheal tubes
- cupped, soft tip
- latex-free
- sterile



FLEXISLIP INTUBATION STYLET

REF.	ORDER SIZE	LENGTH CURVED	O.D.	FOR I.D.	QTY
503700	6	220 mm	1.95 mm	2.5 mm	20
	10	335 mm	3.28 mm	3.5 mm	
	12	366 mm	3.93 mm	4.0 mm	
	14	366 mm	4.60 mm	5.0 mm	

INTUBATION STYLET MADE OF PVC

- for tracheal tubes
- flexible
- cupped, soft tip
- latex-free



INTUBATION STYLET MADE OF PVC

RÜSCH

REF.	ORDER SIZE	LENGTH APPROX.	O.D.	FOR I.D.	QTY
503000	15	350 mm	1.5 mm	2.0–2.5 mm	10
	20	400 mm	2.0 mm	2.5–3.0 mm	
	26	400 mm	2.6 mm	3.0–3.5 mm	
	33	450 mm	3.3 mm	3.5–4.5 mm	
	43	450 mm	4.3 mm	4.5–6.0 mm	
	56	500 mm	5.6 mm	ab 6.0 mm	

- for non-kinking (Oxford) tube
- right-angled
- hard
- latex-free



INTUBATION STYLET MADE OF PVC

RÜSCH

REF.	ORDER SIZE	LENGTH APPROX.	O.D.	FOR I.D.	QTY
502900	1	280 mm	2.6 mm	3.5–5.0 mm	5
	2	370 mm	4.3 mm	5.5–6.0 mm	
	3	410 mm	5.6 mm	7.0–12.0 mm	

ENDOGUIDE TUBE EXCHANGER AND INSERTION AID IN ONE

permanently available large lumen for constant oxygenation of the patient; connection system for 15 mm and Luer-lock connectors, depth markers

- latex-free
- sterile



ENDOGUIDE-T

RÜSCH

REF.	ORDER SIZE	I.D.	O.D.	LENGTH*	QTY
503100	2.5	1.4 mm	2.6 mm	700 mm	1
	6.0	3.2 mm	5.0 mm	830 mm	
503110	6.0	3.2 mm	5.0 mm	525 mm	

* Length without connector

CONNECTOR SET

- Luer-lock connector
- silicone cone
- 15 mm standard connector

