

OROPHARYNGEAL/ NASOPHARYNGEAL AIRWAYS

For quick and safe airway patency without intubation the user can choose from an extensive range of RÜSCH oropharyngeal and nasopharyngeal airways.

All versions have been specifically developed to meet the customer's requirements and thus represent a safe and adequate solution for most indications. RÜSCH offers not only

standard versions made from a variety of materials but also special models for endoscopy that offer simultaneous fixation of the tracheal tube. Thanks to the broad selection, these RÜSCH products have become indispensable in many situations.

OROPHARYNGEAL AIRWAYS

OROPHARYNGEAL AIRWAY		RÜSCH				
DESCRIPTION	REF.	ORDER SIZE	ISO SIZE	COLOUR CODE	QTY	
 <ul style="list-style-type: none"> made of PVC Guedel airway colour-coded clear latex-free 	124700	000	3	colourless	10	
		00	4	pink		
		0	5	blue		
		1	6	black		
		2	7	white		
		3	8	green		
		4	9	yellow		
		5	10	red		
 <ul style="list-style-type: none"> made of polyethylene Guedel airway colour-coded latex-free 	124900 sterile	000	3	colourless	10	
		00	4	pink		
		0	5	blue		
		1	6	black		
		2	7	white		
		3	8	green		
		4	9	yellow		
		5	10	red		
	672901 non-sterile	00	4	pink	50	
		0	5	blue		
		1	6	black		
		2	7	white		
		3	8	green		

OROPHARYNGEAL AIRWAYS

RÜSCH

	DESCRIPTION	REF.	ORDER SIZE	ISO SIZE		QTY
	made of soft rubber for fiberoptic intubation • Guedel airway • slit along upper surface	124400*	3	8		2
			4	9		
			5	10		
 Ref. 124500 Ref. 124501	made of soft rubber • Guedel airway	 124500* black 124501* red	000	3		10
			00	4		
			0	5		
			1	6		
			2	6.5		
			3	8		
			4	9		
			5	10		
			6	12		
	DESCRIPTION	REF.	ORDER SIZE	I.D.	LENGTH APPROX.	QTY
	OPTOSAFE bite block made of PVC, for fiberoptic intubation • Guedel airway • circular inner diameter • latex-free • sterile	 124200	7	7 mm	67 mm	2
			11	11 mm	80 mm	
			13	13 mm	80 mm	
			15	15 mm	95 mm	

*** Attention**
 These products contain natural rubber latex which may cause allergic reactions.

CATH-GUIDE® AIRWAY

HUDSON RCI

	DESCRIPTION	REF.	ORDER SIZE	I.D.	QTY	
	• flexible vinyl with a rigid bite block • Guedel style with three internal channels • available in eight sizes • individually packed	 	41165	5	120 mm	48
			41166	4	110 mm	
			41167	3	100 mm	
			41164	2	90 mm	
			41168	1	80 mm	
			41169	0	70 mm	
			41170	00	60 mm	
			41171	000	55 mm	

BITEGARD® ORAL BITE BLOCK

HUDSON RCI

	DESCRIPTION	REF.	QTY
	• oral bite block  	41140	2