

# SHERIDAN BRONCHIAL TUBES

## SHER-I-BRONCH®

### LEFT-SIDED DOUBLE LUMEN ENDOBRONCHIAL TUBE FOR INTUBATION OF LEFT MAINSTEM BRONCHUS

colour coding distinguishes bronchial from tracheal airways; radiopaque cuff attachment points and tube tip, full-scale depth marks, blue X-ray line for ease of visualisation during fiberoptic scope insertion, tear-resistant tracheal cuff

- sterile
- individually packed



#### SHER-I-BRONCH®

**SHERIDAN**

REF.	ORDER SIZE	2 X I.D.	QTY
5-16028	28	3.78	1
5-16035	35	4.90	
5-16037	37	5.16	
5-16039	39	5.51	
5-16041	41	5.56	

### RIGHT-SIDED DOUBLE LUMEN ENDOBRONCHIAL TUBE FOR INTUBATION OF RIGHT MAINSTEM BRONCHUS

1 large eye (for ventilating right upper lobe bronchial orifice)

unique double cuff design increases the margin of safety, when positioning the tube in the right mainstem bronchus; full-scale depth marks, blue X-ray line for ease of visualisation during fiberoptic scope insertion, tear-resistant tracheal cuff

- sterile
- individually packed



#### SHER-I-BRONCH®

**SHERIDAN**

REF.	ORDER SIZE	2 X I.D.	QTY
5-16128	28	3.78	1
5-16135	35	4.90	
5-16137	37	5.16	
5-16139	39	5.51	
5-16141	41	5.56	

\* packaged separately with every SHER-I-BRONCH® endobronchial tube