

# SHERIDAN ENDOTRACHEAL TUBES

Supplied under the SHERIDAN brand name, our endotracheal tubes offer all the features needed by today's health-care professionals, including a high level of tactile feed-back from the pilot balloon, thermosensitive PVC,

which allows the tubes to soften at body temperature, as well as a smooth tube tip that is moulded and gently cupped to minimise trauma during intubation.

## SHERIDAN ISIS HVT

**ENDOTRACHEAL TUBE WITH INTEGRATED PORT FOR OPTIONAL SUBGLOTTIC SECRETION SUCTIONING WITH HIGH VOLUME TAPERED CUFF**  
nasal/oral, 1 eye (Murphy)

The ideal tube to provide the best of care to all your patients, both over the short- & long-term.

The unique convertible nature of the Teleflex ISIS HVT frees clinicians from the burden of choosing which tube is best for the patient at the time of intubation. ISIS features an integrated suction port and separate suction line allowing for subglottic secretion removal whenever necessary. Both connection ports can be sealed upon disconnection, reducing the risk of cross-contamination when not in use. When needed, the suction tube attaches to the ISIS HVT via a secure locking connection.

- sterile



### SHERIDAN ISIS HVT

**SHERIDAN**

REF.	ORDER SIZE/I.D.	QTY
5-13012	6.0 mm	10
5-13013	6.5 mm	
5-13014	7.0 mm	
5-13015	7.5 mm	
5-13016	8.0 mm	
5-13017	8.5 mm	
5-13018	9.0 mm	



### SUCTION LINE

REF.	QTY
5-23000	5

# SHERIDAN/HVT® TRACHEAL TUBE



SHERIDAN/HVT TRACHEAL TUBE **SHERIDAN**

REF.	ORDER SIZE/I.D.	QTY
5-10310	5.0 mm	10
5-10311	5.5 mm	
5-10312	6.0 mm	
5-10313	6.5 mm	
5-10314	7.0 mm	
5-10315	7.5 mm	
5-10316	8.0 mm	
5-10317	8.5 mm	
5-10318	9.0 mm	
5-10320	10.0 mm	

**HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF**  
nasal/oral, 1 eye (Murphy)

The high-volume, large diameter cuff provides a positive tracheal wallseal, while the low-pressure cuff minimises capillary restriction. The tapered cuff reduces the risk of micro-aspiration, increases ease of use & minimises the risk of herniation.

- sterile

# SHERIDAN/CF® TRACHEAL TUBE



SHERIDAN/CF TRACHEAL TUBE **SHERIDAN**

REF.	ORDER SIZE/I.D.	REF.	ORDER SIZE/I.D.	QTY
CF Magill Type		CF Murphy Eye		10
5-10206	3.0 mm	5-10106	3.0 mm	
	-	5-10107	3.5 mm	
5-10208	4.0 mm	5-10108	4.0 mm	
	-	5-10109	4.5 mm	
5-10210	5.0 mm	5-10110	5.0 mm	
	-	5-10111	5.5 mm	
5-10212	6.0 mm	5-10112	6.0 mm	
5-10213	6.5 mm	5-10113	6.5 mm	
5-10214	7.0 mm	5-10114	7.0 mm	
5-10215	7.5 mm	5-10115	7.5 mm	
5-10216	8.0 mm	5-10116	8.0 mm	
5-10217	8.5 mm	5-10117	8.5 mm	
5-10218	9.0 mm	5-10118	9.0 mm	
5-10219	9.5 mm	5-10119	9.5 mm	
5-10220	10.0 mm	5-10120	10.0 mm	

**CLOSE FITTING, LOW-PRESSURE CUFF**  
nasal/oral, 1 eye (Murphy or Magill)

By providing ease of intubation and increasing visualisation, the opaque, close fitting, low- pressure cuff is ideal for short-term surgical intubations.

- sterile

# SHERIDAN UNCUFFED™ TRACHEAL TUBE



## SHERIDAN UNCUFFED TRACHEAL TUBE

**SHERIDAN**

REF.	ORDER SIZE/I.D.	QTY
5-10404	2.0 mm	10
5-10405	2.5 mm	
5-10406	3.0 mm	
5-10407	3.5 mm	
5-10408	4.0 mm	
5-10409	4.5 mm	
5-10410	5.0 mm	
5-10411	5.5 mm	
5-10412	6.0 mm	
5-10413	6.5 mm	
5-10414	7.0 mm	

### MADE OF PVC, WITHOUT CUFF

nasal/oral, 1 eye (Murphy)  
cut on (\*) for oral use

Special depth lines at the distal tip indicate the depth of intubation below the chords. The smooth, moulded, gently cupped tip minimises trauma during intubation, while the tube softens at body temperature to conform to the airway anatomy.

- sterile

# SHERIDAN PREFORMED™ TRACHEAL TUBE

## HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF

nasal/oral, 1 eye (Murphy) with cuffed version  
 nasal/oral, 2 eyes (Murphy) with uncuffed version

Sheridan Preformed Tracheal Tubes for oral and maxillofacial surgery allow the anaesthesia circuit connection to be positioned out of the surgical field.

Features include: bold marks at the center of the bend with distance to distal tip indicated, high-volume, large diameter cuff for positive tracheal wall seal, low-pressure cuff to minimise capillary restrictions.

Uncuffed tracheal tubes have two Murphy eyes for enhanced patient safety.

- sterile



### SHERIDAN PREFORMED TRACHEAL TUBE, ORAL



REF.	ORDER SIZE/I.D.	REF.	ORDER SIZE/I.D.	QTY
Uncuffed oral		Cuffed oral		10
5-22006	3.0 mm		-	
5-22007	3.5 mm		-	
5-22008	4.0 mm	5-22208	4.0 mm	
5-22009	4.5 mm	5-22209	4.5 mm	
5-22010	5.0 mm	5-22210	5.0 mm	
5-22011	5.5 mm	5-22211	5.5 mm	
5-22012	6.0 mm	5-22212	6.0 mm	
5-22013	6.5 mm	5-22213	6.5 mm	
5-22014	7.0 mm	5-22214	7.0 mm	
		5-22215	7.5 mm	
		5-22216	8.0 mm	
		5-22217	8.5 mm	
		5-22218	9.0 mm	



### SHERIDAN PREFORMED TRACHEAL TUBE, NASAL



REF.	ORDER SIZE/I.D.	REF.	ORDER SIZE/I.D.	QTY
Uncuffed nasal		Cuffed nasal		10
5-22106	3.0 mm		-	
5-22107	3.5 mm		-	
5-22108	4.0 mm		-	
5-22109	4.5 mm		-	
5-22110	5.0 mm		-	
5-22111	5.5 mm		-	
5-22112	6.0 mm	5-22312	6.0 mm	
5-22113	6.5 mm	5-22313	6.5 mm	
5-22114	7.0 mm	5-22314	7.0 mm	
		5-22315	7.5 mm	
		5-22316	8.0 mm	

# SPIRAL-FLEX REINFORCED TRACHEAL TUBE

## HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF

nasal/oral, with cuffed version

nasal/oral, 1 eye (Murphy) with uncuffed version

reduced risk of kinking due to stainless steel, spiral wound reinforcing wire within the tube wall; 15 mm connector is permanently bonded within the tube, blue cuff

- sterile

### CUFFED ORAL

cuffed, oral precut length with reinforcing sleeve to minimise the risk of patient biting and occluding the tube; the oral length tube has a reinserted SHER-I-SLIP stylet.



### SPIRAL-FLEX REINFORCED TRACHEAL TUBE

**SHERIDAN**

REF.	ORDER SIZE/I.D.	REF.	ORDER SIZE/I.D.	REF.	ORDER SIZE/I.D.	QTY
Cuffed oral		Cuffed oral/nasal		Uncuffed oral/nasal		5
-	-	-	-	5-12906	3.0 mm	
-	-	-	-	5-12907	3.5 mm	
-	-	-	-	5-12908	4.0 mm	
-	-	-	-	5-12909	4.5 mm	
5-12510	5.0 mm	5-12610	5.0 mm	5-12910	5.0 mm	
5-12511	5.5 mm	5-12611	5.5 mm	-	-	
5-12512	6.0 mm	5-12612	6.0 mm	-	-	
5-12513	6.5 mm	5-12613	6.5 mm	-	-	
5-12514	7.0 mm	5-12614	7.0 mm	-	-	
5-12515	7.5 mm	5-12615	7.5 mm	-	-	
5-12516	8.0 mm	5-12616	8.0 mm	-	-	
5-12517	8.5 mm	5-12617	8.5 mm	-	-	
5-12518	9.0 mm	5-12618	9.0 mm	-	-	

**SHERIDAN LTS™**  
**TRACHEAL TUBE FOR MICROLARYNGEAL**  
**TRACHEAL SURGERY**  
 nasal/oral, 1 eye (Murphy)

small O.D. tube allows greater access to the surgical field, yet provides sufficient ventilation of the patient; same tube length and cuff diameter as an 8.0mm cuffed tracheal tube; high-volume, yellow cuff is more visible and provides a low-pressure tracheal wall seal

- sterile
- individually packed



SHERIDAN LTS			<i>SHERIDAN</i>
REF.	ORDER SIZE/I.D.	SIZE/O.D.	QTY
5-11108	4.0 mm	5.8 mm	10
5-11110	5.0 mm	7.1 mm	
5-11112	6.0 mm	8.5 mm	