


SHER-I-SLIP® INTUBATION STYLET

dependable stiffness retains desired curve during oral intubation; blue colour improves visualisation; special material reduces stylet-to-tube friction; smooth moulded distal tip end



SHER-I-SLIP INTUBATION STYLET *HUDSON RCI*

REF.	LENGTH APPROX.	O.D.	FOR TUBE SIZES I.D.	QTY
5-15101	340 mm	10 Fr	5.0 - 10.0 mm	25
5-15102	250 mm	5 Fr	2.0 - 4.5 mm	
5-15103	X-Long	10 Fr	5.0 - 10.0 mm	

SOFT-TIP™ INTUBATION STYLET

atraumatic tip to minimise trauma; slippery outer finish makes for easy transition into endotracheal tube



SOFT-TIP INTUBATION STYLET *HUDSON RCI*

REF.	LENGTH APPROX.	O.D.	FOR TUBE SIZES I.D.	QTY
5-15120	270 mm	6 Fr	2.5 - 4.5 mm	20
5-15110	330 mm	10 Fr	4.5 - 6.5 mm	
5-15100	310 mm	14 Fr	6.0 - 10.0 mm	