

# NASOPHARYNGEAL AIRWAYS

## NASOPHARYNGEAL AIRWAY, WENDL PATTERN, MADE OF WIRUPREN®

- with adjustable flange
- latex-free
- sterile



## NASOPHARYNGEAL AIRWAY, WENDL PATTERN, MADE OF SOFT RUBBER

- with adjustable flange








**Attention:**  
These products contain natural rubber latex which may cause allergic reactions.

## NASOPHARYNGEAL AIRWAY MADE OF PVC


- transparent
- sterile



NASOPHARYNGEAL AIRWAYS					RÜSCH
REF.	ORDER SIZE	I.D.	O.D.	LENGTH	QTY
185420  	12	2,0 mm	4.00 mm	95 mm	5
	14	2,5 mm	4.67 mm	95 mm	
	16	3,0 mm	5.33 mm	95 mm	
	18	3,5 mm	6.00 mm	95 mm	
	20	4,5 mm	6.67 mm	125 mm	
	22	5,0 mm	7.33 mm	170 mm	
	24	5,5 mm	8.00 mm	170 mm	
	26	6,0 mm	8.67 mm	170 mm	
	28	6,5 mm	9.33 mm	170 mm	
	30	7,0 mm	10.00 mm	170 mm	
	32	7,5 mm	10.67 mm	170 mm	
	34	8,0 mm	11.33 mm	170 mm	
	36	8,5 mm	12.00 mm	170 mm	
125200 	12	2.0 mm	4.0 mm	95 mm	5
	14	2.5 mm	4.7 mm	95 mm	
	16	3.0 mm	5.3 mm	95 mm	
	18	3.5 mm	6.0 mm	95 mm	
	20	4.0 mm	6.7 mm	115 mm	
	22	4.5 mm	7.3 mm	125 mm	
	24	5.0 mm	8.0 mm	170 mm	
	26	5.5 mm	8.7 mm	170 mm	
	28	6.0 mm	9.3 mm	170 mm	
	30	6.5 mm	10.0 mm	170 mm	
	32	7.0 mm	10.7 mm	170 mm	
	34	7.5 mm	11.3 mm	170 mm	
	36	8.0 mm	12.0 mm	170 mm	
125410  	12	2.5 mm	3.97 mm	60 mm	10
	14	3.0 mm	4.64 mm	72 mm	
	16	3.5 mm	5.30 mm	85 mm	
	18	4.0 mm	5.97 mm	105 mm	
	19	4.5 mm	6.29 mm	110 mm	
	20	5.0 mm	6.63 mm	115 mm	
	22	5.5 mm	7.29 mm	125 mm	
	24	6.0 mm	7.96 mm	130 mm	
	26	6.5 mm	8.63 mm	140 mm	
	28	7.0 mm	9.29 mm	155 mm	
	30	7.5 mm	9.96 mm	165 mm	
	32	8.0 mm	10.63 mm	170 mm	
	34	8.5 mm	11.29 mm	175 mm	
36	9.0 mm	11.95 mm	180 mm		

**NASOPHARYNGEAL AIRWAY,  
MADE OF SILKOLATEX®  
RÜSCH GOLD**

- with adjustable flange and widened end
- sterile




**NASOPHARYNGEAL AIRWAY,  
ROBERTAZZI PATTERN,  
MADE OF WIRUPREN®**

- latex-free
- sterile



**NASOPHARYNGEAL AIRWAY,  
WENDL PATTERN, MADE OF  
SOFT RUBBER**

- double lumen
- with adjustable flange
- oxygen insufflation tube



**NASOPHARYNGEAL AIRWAYS**

**RÜSCH**

REF.	ORDER SIZE	I.D.	O.D.	LENGTH	QTY
185200* ⊗	20	4.5 mm	6.7 mm	125 mm	10
	22	5.0 mm	7.3 mm	170 mm	
	24	5.5 mm	8.0 mm	170 mm	
	26	6.0 mm	8.7 mm	170 mm	
	28	6.5 mm	9.3 mm	170 mm	
	30	7.0 mm	10.0 mm	170 mm	
	32	7.5 mm	10.7 mm	170 mm	
	34	8.0 mm	11.3 mm	170 mm	
	36	8.5 mm	12.0 mm	170 mm	
	185410 ⊗	20	3.6 mm	6.67 mm	
22		4.6 mm	7.33 mm	105 mm	
24		5.2 mm	8.00 mm	105 mm	
26		5.8 mm	8.67 mm	110 mm	
28		6.5 mm	9.33 mm	117 mm	
30		7.0 mm	10.00 mm	118 mm	
32		7.6 mm	10.67 mm	140 mm	
34		7.9 mm	11.33 mm	155 mm	
36		9.0 mm	12.00 mm	170 mm	
125600* ⊗		24	5.5 mm	8.0 mm	160 mm
	26	6.0 mm	8.7 mm	160 mm	
	28	6.5 mm	9.3 mm	170 mm	
	30	7.0 mm	10.0 mm	170 mm	
	32	7.6 mm	10.7 mm	170 mm	

**\*Attention**

*These products contain natural rubber latex which may cause allergic reactions.*



**STANDARD CONNECTOR**

**RÜSCH**

REF.	ORDER SIZE/I.D.	DESCRIPTION	QTY
501003	1.5 mm – 11.0 mm	made of plastic, latex-free, 15 mm connection (O.D.)	10