



RÜSCH RECTAL/PHARYNGEAL TEMPERATURE PROBE

For reliable and easy monitoring of the core body temperature




AT A GLANCE

- measuring accuracy of $\pm 0,2^{\circ}\text{C}$
- measurement range $25^{\circ}\text{C} - 45^{\circ}\text{C}$
- graduated
- smooth soft material
- easy insertion
- latex-free
- DEHP-free
- clean
- single use

OPTIMUM PERFORMANCE AND CONVENIENCE

The rÜSCH rectal-/pharyngeal temperature probe provided by Teleflex is a temperature probe and indicated for measuring the core body temperature during routine temperature management on patients in OR (operating room) and ICU (intensive care unit).

The temperature sensor of this device has a measuring accuracy of $\pm 0,2^{\circ}\text{C}$ over a range of $25^{\circ}\text{C} - 45^{\circ}\text{C}$. The probe is made of extraordinarily soft and smooth material and can thus be easily inserted. The material does not contain any latex nor DEHP.








- ① connector with ergonomic grip for safe and simple handling
- ② graduated
- ③ catheter shaft made of extremely soft and smooth material
- ④ connecting cable

RECTAL/PHARYNGEAL TEMPERATURE PROBE

REF.	LENGTH	QTY
1016	<ul style="list-style-type: none"> • introduction length catheter: 46 cm • total length catheter: 73 cm 	30

THE ADAPTER CABLES – MALE CONNECTION

REF.	DESCRIPTION	LENGTH	QTY
The adapter cables are compatible with all monitors meeting the type series 400 industrial standard.			
	81-101400 Hospital plug, straight, for Datex Cardiocap, Datex-Ohmeda Astiva, Datex-Ohmeda Anesthesia, as well as for GE Hellige SMK, GE Eagle Dash, GE Solar 8000 M, GE 8000 Solo, GE 7000, GE Tram, Nihon Kohden Life Scope, as well as for Datascope Passport, Datascope Point of View, Datascope 200/2000A, 800, and Spacelaps 90600, 90400, 90300	5 m	1
	81-1010SM Siemens round plug, 7-pole, for Dräger (Siemens): Infinity Delta XL, SC 6000, SC 7000, SC 8000, SC 9000	5 m	1
	81-1010HP HP round plug, 2-pole, for Phillips Intellie Vue, Sure Signs - VM Serie, MP 2, MP 5, MP 20, MP 30, MP 40, MP 50, MP 60, MP 90 as well as for Hewlett Packard Viridia 24 C, 78300, Merlin	5 m	1

DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business & Technology Park · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

WWW.TELEFLEX.COM

94 06 22 - 00 00 01 · REV A · 10 11 02 · All data current at time of printing. Subject to technical changes without further notice.

teleflex