

DISPOSABLE LARGE VOLUME NEBULISER FOR MECHANICAL VENTILATION

- 500 ml
- defined maximum and minimum water level lines
- concentration: 28%-98%
- latex-free
- single use
- individually packed

MATERIALS

- housing: polyethylene
- adaptor: polycarbonate



DISPOSABLE LARGE VOLUME NEBULISER FOR MECHANICAL VENTILATION		<i>HUDSON RCI</i>
REF.		QTY
41770		50

AQUATHERM® III EXTERNAL ADJUSTABLE ELECTRONIC HEATER

- attaches onto the nebuliser system to provide heated aerosol
- adjustable aerosol temperature to meet patient needs with a sensitive temperature controller
- robust, durable heater

AQUATHERM III EXTERNAL ADJUSTABLE ELECTRONIC HEATER		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
050-12 D	220 - 240 V, 50 Hz	1
050-12 E	United Kingdom style power plug	

