

TRACHEOSTOMY MASK

- for tracheostomy and laryngectomy aerosol therapy
- tubing connector swivels 360° and accepts 22 mm I.D. corrugated tubing
- clean
- single use
- individually packed



TRACHEOSTOMY MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
1075	with adult size mask	50
41076	with paediatric size mask	

MATERIALS

- mask: PV with metal bracket
- adaptor: polyethylene
- neck band: polychloroprene

AEROSOL MASK

- clear, soft vinyl
- specifically designed for aerosol therapy
- connector accepts 22 mm aerosol tubing
- clean
- single use
- individually packed



ELONGATED AEROSOL MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
1083	with adult size mask	50
41085	with paediatric size mask	

FACE TENT

- clear, soft vinyl
- for high-humidity aerosol therapy
- connector accepts 22 mm aerosol tubing
- clean
- single use
- individually packed



FACE TENT

HUDSON RCI

REF.	DESCRIPTION	QTY
41095		50