

# INCENTIVE SPIROMETERS

## THE EFFECTIVE BREATHING THERAPY WITH SUSTAINED MAXIMAL INSPIRATION

The incentive spirometers from HUDSON visualise patients' breathing performance (ball, markings) and thus induce them to take deeper and longer breaths. Since they do not need to be reprocessed, these single-patient products are also very economical.

## TRIFLO® II

- flow-oriented, 3-ball incentive spirometer with mouthpiece and tubing
- three colour-coded balls/three chambers
- minimum flow imprinted on each chamber
- compact design and break-resistant plastic
- latex-free

TRIFLO II		HUDSON RCI	
REF.	FLOW RATE	SIZE	QTY
8884717395	600-1200 ml/sec.	housing (length/width/height) 13.5/7/14 cm tubing 28 cm long mouthpiece 6.5 cm long	12



## VOLDYNE®

- volume-oriented incentive spirometer with mouthpiece and tubing
- two sizes to cater for small as well as large inspired volumes
- compact, economical design with built-in handle
- markings on both sides
- latex-free

VOLDYNE 2500/VOLDYNE 5000		HUDSON RCI		
REF.	DESCRIPTION	VOLUME MEASUREMENT	SIZE	QTY
8884719041	Voldyne 2500 for adult and paediatric use	250 ml to 2500 ml	housing (length/width/height) 13.5/7/14 cm tubing 28 cm long mouthpiece 6.5 cm long	12
8884719033	Voldyne 5000 for adult use	250 ml to 5000 ml	housing (length/width/height) 13.5/7/14 cm tubing 28 cm long mouthpiece 6.5 cm long	



Ref. 8884719041

## LUNG VOLUME EXERCISER

- flow-oriented incentive spirometer with mouthpiece and tubing
- dual-chamber design creates a constant resistance for deep and prolonged inspiration
- break-resistant material
- latex-free

LUNG VOLUME EXERCISER  
FOR ADULT AND PAEDIATRIC USE

*HUDSON RCI*

REF.	FLOW RATE	QTY
41750	200-1200 ml/sec. (adjustable flow settings)	10

