

AQUAPAK®

PREFILLED NEBULISERS


AQUAPAK prefilled nebulisers offer versatile, sterile systems to provide cold or heated nebulisation. Choose from a wide variety of reservoirs, adaptors and solutions to meet specific patient requirements efficiently and economically.

- latex-free
- gamma sterile




- select between two AQUAPAK precision nebuliser adaptors for delivery of a consistent aerosol over a wide range of oxygen concentrations
- AQUATHERM III heater provides heated aerosol treatment
- prefilled AQUAPAK sterile reservoirs are available in four convenient sizes: 440 ml, 760 ml, 1070 ml, and 2200 ml
- available with both reservoir and sterile adaptor conveniently packed together
- AQUAPAK reservoir may be used with other HUDSON humidification and nebulisation products, making the system truly versatile

NEBULISER ADAPTOR 028 HUDSON RCI

| REF. | DESCRIPTION | QTY |
|--------|--|-----|
| 403128 |  <ul style="list-style-type: none"> • Venturi-style entrainment ring allows the selection of oxygen concentrations from 28% to 98% • adaptor: polycarbonate • individually packed • sterile | 50 |

NEBULISER ADAPTOR 033 HUDSON RCI

| REF. | DESCRIPTION | QTY |
|--------|--|-----|
| 403133 |  <ul style="list-style-type: none"> • quiet adaptor allows the selection of oxygen concentrations from 33% to 98% to meet the needs of the ICU patient • adaptor: polycarbonate • individually packed • sterile | 50 |

CLOSED STERILE WATER SYSTEM FOR MECHANICAL NEBULISATION + ADAPTOR HUDSON RCI

| REF. | 440 ML | 760 ML | 1070 ML | 2200 ML | STERILE WATER | 0.45% SALINE | 0.9% SALINE | ADAPTOR 028 | ADAPTOR 033 | QTY |
|--------|--------|--------|---------|---------|---------------|--------------|-------------|-------------|-------------|-----|
| 400400 | x | | | | x | | | | | 20 |
| 403128 | | | | | | | | x | | 50 |
| 403133 | | | | | | | | | x | 50 |
| 403700 | | x | | | x | | | | | 10 |
| 403705 | | x | | | | x | | | | 10 |
| 403709 | | x | | | | | x | | | 10 |
| 403728 | | x | | | x | | | x | | 10 |
| 404000 | | | x | | x | | | | | 10 |
| 404128 | | | x | | x | | | x | | 10 |
| 404200 | | | | x | x | | | | | 4 |
| 404290 | | | | x | | | x | | | 4 |
| 404428 | x | | | | x | | | x | | 20 |
| 404739 | | | | x | | | x | x | | 4 |