

# AQUAPAK®

## PREFILLED HUMIDIFIERS

HUDSON AQUAPAK® prefilled humidifiers combine an easy-to-use adaptor with a choice of two sterile water reservoir bottles. The humidifier adaptor incorporates an audible alarm that alerts the clinician to flow restriction or occlusion of the humidifier or oxygen tubing. The reservoir bottles have a built-in oxygen tubing connector with an easy snap-off trigger design. A micro-diffuser, moulded into the reservoir bottle, produces smaller bubbles and increases surface area agitation allowing for noise-free humidification. HUDSON AQUAPAK prefilled humidifiers provide performance, safety and value. The complete set with its insufflation adaptor facilitates handling.

- latex-free
- gamma sterile
- single use
- individually packed

**MATERIALS**

- bottle: polypropylene
- adaptor: polycarbonate
- tubing: PVC
- valve: silicone



### AQUAPAK STERILE WATER, HUDSON RCI INCLUDES HUMIDIFIER ADAPTOR

REF.	DESCRIPTION	QTY
400340	340 ml sterile aqua dest., includes Ref. 400040 sterile Humidifier Adaptor	20
400640	650 ml sterile aqua dest., includes Ref. 400040 sterile Humidifier Adaptor	10



### AQUAPAK® STERILE WATER, HUDSON RCI WITHOUT HUMIDIFIER ADAPTOR

REF.	DESCRIPTION	QTY
400301	340 ml sterile aqua dest., water reservoir only, for use with Ref. 400040 Humidifier Adaptor	20
400601	650 ml sterile aqua dest., water reservoir only, for use with Ref. 400040 Humidifier Adaptor	10

### AQUAPAK STERILE WATER HUDSON RCI

REF.	DESCRIPTION	QTY
400400	440 ml sterile aqua dest., water reservoir only, for use with Ref. 400011 Silent Adaptor	20
403700	760 ml sterile aqua dest., water reservoir only, for use with Ref. 400011 Silent Adaptor	10

### AQUAPAK SILENT ADAPTOR HUDSON RCI

REF.	DESCRIPTION	QTY
400011	<ul style="list-style-type: none"> <li>• no bubbles, vibration or noise</li> <li>• pressure relief valve, 140-210 cm H<sub>2</sub>O</li> <li>• audible pop-off</li> <li>• designed for use with Ref. 400400 or 403700 sterile water</li> <li>• optional external heater capability for use with Ref. 050-12 AQUATHERM series heaters</li> <li>• single use</li> </ul>	50



### HUMIDIFIER ADAPTOR HUDSON RCI

REF.	DESCRIPTION	QTY
400040	<ul style="list-style-type: none"> <li>• sterile, clear adaptor body for use with Ref. 400301 or Ref. 400601</li> <li>• pressure relief valve 350-700 cm H<sub>2</sub>O</li> <li>• audible pop-off</li> <li>• single use</li> </ul>	120

