

NASAL CANNULAS

PATENTED NASAL CANNULAS WITH OPTIMAL PERFORMANCE

The curved nasal tips are thinner and softer to the touch and follow the anatomical contours of the patient’s nose. The over-the-ear style cannulas are available in a variety of configurations offering the customer a choice of nasal tips, oxygen tubing, and tubing connectors.

- flow rate: 1-7 LPM
- latex-free
- clean
- single use



Ref. 1103



Ref. 1104

ADULT OVER-THE-EAR NASAL CANNULAS

HUDSON RCI

REF.	NASAL TIP		TUBING			STANDARD TUBING CONNECTOR	QTY
	FLARED	NON-FLARED	STANDARD	STAR-LUMEN	LENGTH		
1103		x	x		210 cm	x	50
1104	x		x		210 cm	x	
41110		x		x	210 cm	x	

SOFTECH® NASAL CANNULAS

Extremely lightweight and flexible, the SOFTECH cannulas are anatomically designed to enhance patient acceptance and comfort. The over-the-ear style optimises fit and stability. Paediatric and infant sizes are also available.

- latex-free
- clean
- single use

HUDSON cannulas are available with either standard or universal hose end connectors. The universal connectors are convenient and economical because they attach directly to flow meters, thus eliminating the need for nipple-and-nut adaptors.



Ref. 41820


SOFTECH NASAL CANNULAS

HUDSON RCI

REF.	SIZE	TUBING		TUBING CONNECTOR		SUGGESTED FLOW RATE	QTY
		STAR-LUMEN	LENGTH	STANDARD	UNIVERSAL		
41820	Adult	x	210 cm	x		4 LPM	50
41826	Paediatric	x	210 cm		x	3.5 LPM	
41828	Infant	x	210 cm	x		1.5-2 LPM	

MATERIAL
• PVC

OXYGEN SUPPLY TUBING



HUDSON offers both standard (fig. 1) and kink- and crush-resistant STAR LUMEN (fig. 2) tubing. The distinctively designed tubing features a 5-channel inner lumen for continuous oxygen flow.

HUDSON oxygen supply tubing is available with either standard (fig. 3) or universal hose end connectors (fig. 4). The universal connectors are convenient and economical because they attach directly to flow meters, thus eliminating the need for nipple-and-nut adaptors.

- latex-free
- clean
- single use

MATERIAL
• PVC

OXYGEN SUPPLY TUBING

HUDSON RCI

REF.	TUBING		CONNECTOR		LENGTH	QTY
	STANDARD	STAR-LUMEN	STANDARD	UNIVERSAL		
41113		x	x		210 cm	50
1115	x		x		210 cm	50
41118		x	x		420 cm	50
41119		x	x		760 cm	25
41120		x	x		1520 cm	25
41925	x			x	210 cm	50