

DISPOSABLE HUMIDIFIER

HUDSON offers a choice of humidifiers that can be used with the clinician’s water source. Available in disposable form, these humidifiers incorporate an audible pressure-relief valve to signal occlusion of the humidifier or oxygen tubing.

- disposable, durable plastic with preset audible alarm at 280 cm H₂O (4 psi) pressure relief
- defined minimum and maximum water level lines
- latex-free
- single use
- individually packed



MATERIALS

- housing: polyethelene
- adaptor: polypropylene
- pressure compensation valve: silicone

DISPOSABLE HUMIDIFIER

HUDSON RCI

REF.	FILLING VOLUME	QTY
3230	500 ml	50